



SUITABLE FOR LABORATORY AUTOMATION

**BC+™**

**POSITIVE HEMOCULTURE  
PROCESSING DEVICE**



Innovating Together™

[copangroup.com](http://copangroup.com)



# BC+™

BC™ + is a vacuum tube designed to transfer and transport positive hemoculture for subsequent culture or microscopic analysis

READY TO USE

VISUALLY VERIFIABLE SAMPLE VOLUME ASPIRATION

SAMPLE TRANSPORT CONDITIONS: 8H AT RT

EASY TRANSFER OF POSITIVE HEMOCULTURE


SUITABLE FOR AUTOMATIC PROCESSING (WASP™)



## HOW TO USE

Handle the tube under a laminar flow cabinet

- > **Clean** the hemoculture and the BC+™ rubber stoppers with appropriate disinfectant
- > Using a validated piercer, **pierce the cap** of the hemoculture bottle
- > Obliquely holding the bottle and the piercer in one of your hands, with the other hand **insert** the BC+™ in the piercer allowing the BC+™ vacuum aspirate the hemoculture
- > After hemoculture is **transferred** in BC+™ tube, turn bottle and BC+™ in horizontal position and remove the BC+™ tube from the piercer
- > Carefully **remove** the piercer from the hemoculture tube
- > **Check** that the drawn volume of hemoculture in BC+™ tube fits your analysis requirements
- > **Use** the sample in BC+™ for bacterial culture or microscopic analyses



## ORDERING INFORMATION

Code	Product Description	Pack Size
2U026S	BC+™ vacuum tube 12x80 mm, red cap	600 pieces - 12 resealable bags of 50 pieces

Please contact [customer-care@copangroup.com](mailto:customer-care@copangroup.com) for complete ordering information

This document contains information on products which is targeted to a wide range of audiences and may contain product details or information otherwise not accessible or valid in your country. Please be aware that we do not take any responsibility for accessing such information which may not comply with any valid legal process, regulation, registration or usage in the country of your origin.



COPAN ITALIA SPA

Via F. Perotti 10  
25125 - Brescia,  
Italy

t. +39 030 2687211  
f. +39 030 2687250  
[info@copangroup.com](mailto:info@copangroup.com)  
[www.copangroup.com](http://www.copangroup.com)

Brochure hardcopies may not include the latest updates and changes. Please refer to Copan website to view and/or download the most recent controlled version of the brochure. Brochures are mainly intended for marketing purposes. Always consult product inserts, instructions for use and manuals for the appropriate use of the products. ©2016 Copan group. All rights reserved. The trademarks mentioned herein are property of Copan.